

Automated RNA Extraction Machine MaxExtraction96 (Throughput-96)

Introduction: This automated extraction machine based on a magnetic rod-based nucleic acid separation technology in combination with extraction reagents based on preloaded magnetic beads process 1 to 96 samples and purify nucleic acids from a variety of materials, such as: blood, cells, and viruses. To achieve fully automated nucleic acid purification, magnetic beads are absorbed, transported, and released by a specific magnetic rod.

Features



Strong Magnetic Force

Magnetic beads are recovered at a rate of ≥ 98 percent using a 5500 Gauss magnetic rod.



Customer Support

Our team supports customers creating and editing protocols in addition to pre-programmed protocols to fulfill a variety of needs.



Convenient

It fits into a standard biological safety cabinet.



Accurate Temperature Control

Heating for lysis and elution is done automatically, with a quick heating speed.



Prevents Cross-contamination

The UV sterilizing module and intelligent magnetic rod motion control system efficiently eliminate cross-contamination between wells.



Smart Door

The application automatically pauses when the door is opened and resumes when the door is closed.



Flexible

Suitable for use with a variety of magnetic bead-based extraction reagents




Fast Operating Time

Within few minutes this machine can extract RNA or DNA from up to 96 samples.

Automated RNA Extraction Machine MaxExtraction96 (Throughput-96)

Outstanding Magnetic Rod

The driving device is connected with a high-performance motor and a magnetic rod with a considerable vibration amplitude. In order to provide a good and uniform mixing effect, the vibration amplitude can be adjusted according to the solution volume. The rod is controlled by a ball screw actuator, which ensures smooth operation, high precision, and extended service life. To prevent instrument failure, each moving component is safeguarded by a limit position protection mechanism.




Strong Magnetic Rod


The magnetic beads are adsorbed on the head of the magnetic bar, ensuring that the elution buffer can still cover all of the magnetic beads even with a limited elution volume. The high yield of nucleic acids is ensured by the excellent recovery of magnetic beads.

Strong Magnetic Rod

Ordinary Magnetic Rod



Automated RNA Extraction Process



Samples To release the nucleic acid, samples are lysed in Lysis buffer.	Transfer the magnetic beads to the lysis buffer, fully blend, and make sure that the nucleic acids are adsorbed onto the appropriate coated material on the magnetic beads' surface.	Clean the surface of magnetic beads to eliminate any excess protein or salt.	The magnetic beads are transferred to the elution buffer to be mixed thoroughly, the nucleic acid falls off the surface of the magnetic beads and dissolves into the elution buffer.
-----------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Specifications

Product name	MaxExtraction96
Model	A11
Certification	CE/ RoHS
Extraction channel	1-96
Consumables	96 well plate + Tip comb
Nucleic acid extraction time	15-35 minutes
Temperature control precision	0.5°C
Temperature control accuracy	±1.5°C
Temperature uniformity	±1.0°C

Heating range	Room temperature ~95°C
Nucleic acid extraction purity	1.8≤OD260/OD280≤2.0
Inter-well purification variation	CV<3%
Magnetic beads recovery	≥98%
Touch screen size	7 inch color touch screen
Disinfection/decontamination method	UV
Input power	AC 100-240V~, 5.9-2.7A, 50/60Hz
Product weight	52±1kg
Product size	743mm*465mm*447mm

Automated RNA Extraction Machine MaxExtraction96 (Throughput-96)

Viral RNA Extraction Kit – MaxExtraction96 Kit



This Viral RNA extraction kit is compatible with MaxExtraction96 (A11) automated nucleic acid extraction system for extracting ADN/ARN genomically from microorganisms pathogens in samples such as serum, plasma, cultured cells, saliva, alveolar lavado liquid, aspirates, and nasofarangeos.

Name

Viral RNA Extraction Kit by Magnetic Beads

Cat. No.: A12

Size: 96 Reactions

Shipping temperatura: Room temperature

Storage temperatura:

Room temperature.

Fort the long term keep between: 2-8°C

Expiry date: 12 months

Compatibility:

MaxExtraction96 RNA extraction machine (A11)

Advantages

Stable Proteinase K

Our Protease K is stable at ambient temperature, this kit can be transported without ice.

Strong Binding

Magnetic beads designed for pathogenic microbial genomes offer outstanding nucleic acid adsorption properties.

Strong Sealing

High viscosity sealing tape prevents fluid leakage.

Highly Sensitive

This kit can recover trace DNA/RNA from harmful bacteria with good reproducibility and purity.

Procedures



Tear off the sealing film



Add the sample and proteinase K



Load into the instrument



The nucleic acid extraction is completed

Automated RNA Extraction Machine MaxExtraction96 (Throughput-96)

Components

Components	Quantity	Capacity	Components
Tip comb	1	/	/
Plate for samples	1	500 μ L \times 96	Guanidine Hydrochloride, TritonX-100, EDTA, etc
Beads plate	1	200 μ L \times 96	Magnetic Beads
Washing 1 plate	1	600 μ L \times 96	Guanidine Hydrochloride, EDTA, etc
Washing 2 plate	2	600 μ L \times 96	75% ethanol
Elution plate	1	100 μ L \times 96	TE buffer
Proteinase K	1mL x 2	/	Proteinase K

Dimensions

